

INTESTINAL ORGANOID DRUG SCREENING



Purpose

The SOP-ADSI-6.0 was issued to describe the timeline, workflow and procedure of anti-tumor activity testing of drugs on wild-type versus tumor organoids.

Scope

SOP-ADSI-6.0 is intended to cover all resources, personnel and equipment needed to investigate the response (viability) of wild-type as well as tumor organoids on treatment with anti-tumor drugs by using the RealTime-Glo MT Cell Viability Assay.

Introduction

Therapeutic strategies that are based on the individual genetic profile of a patient represent a new frontier of applied cancer research. Especially in cancer research, standard cell culture conditions fail to properly mimic the parental tumor architecture and microenvironment. In this context, tumor-organoids are of special relevance. Tumor-organoids are three-dimensional in vitro cultures derived from the cells obtained from the patient tumor. Tumor-organoids keep the special properties of the original patient's tumor with its microenvironment and mirror the key-features. Thus, tumor-organoids are an ideal tool to identify patient-specific therapies by performing drug-screenings on this primary patient material. The aim of this SOP is to provide protocols for the analysis of tumor-organoids and organoids viability in response to anti-cancer drugs.

1. Equipment, cell culture media, reagents and solutions

1.1 Equipment

- Incubator Binder APT.line™ (150CE2)
- Laminar Air-Flow Labculture Plus ESCO, Class II BSC

- Micro Plate Reader (Mithras LB940)
- Centrifuge (4°C) Hettich Zentrifugen ROTIXA A50 RS
- ImageXpress Micro XLS Widefield High-Content Analysis System
- Table centrifuge VWR Mini Star
- Analog Vortex Mixer (VWR)
- Waterbath

1.2

Materials

- 24 well Culture plates (Greiner bio-one #662160)
- 96-well-plate (Corning #3603)
- Sterile Erlenmeyer flask (200 ml, 500ml)
- Sterile 15 ml and 50 ml tubes
- Sterile 0,5 ml, 1,5 ml, 2 ml and 5 ml tubes
- Set of pipettes 10 µl, 20 µl, 200 µl, 1000 µl and pipette tips
- Pipette Aid and pipettes (10 ml, 25 ml, 50 ml)
- Box with crashed ice

1.3

Chemicals

- Advanced DMEM/F-12 (Gibco-Thermo Fisher #12634028)
- GlutaMax 100 x (GIBCO #35050038)
- HEPES (1 M) (GIBCO #15630080)
- Pen/Strep 100 x (GIBCO #15070063)
- 0.25 % Trypsin-EDTA (Sigma-Aldrich #T4049-100ML)
- WCM* (Broutier 2016)
- RCM** (Broutier 2016)
- NCM*** (Broutier 2016)
- Matrigel (Corning #356231)
- DMSO Hybri-Max (Sigma #D2650)
- Nicotinamide (Sigma-Aldrich #N0636-100G)
- N-Acetylcyteine (Sigma-Aldrich #A9165-25G)
- B-27 (Thermo Fisher #17504044)

* Wnt3a conditioned medium was obtained from L-Wnt3a cells. MTA for the use of producer cell line was obtained from the Hubrecht Institute, Uppsalalaan 8, 3584 CT Utrecht, The Netherlands

** R-spondin1-conditioned medium was produced from 293T-HA-Rspol-Fc producer cell line. MTA for the use of producer cell line was obtained from the Hubrecht Institute, Uppsalalaan 8, 3584 CT Utrecht, The Netherlands

*** Noggin-conditioned medium from HEK293-mNoggin-Fc cell line. MTA for the use of producer cell line was obtained from the Hubrecht Institute, Uppsalalaan 8, 3584 CT Utrecht, The Netherlands

- A83-01 (R&D Systems Europe #2939)
- SB 202190 (Sigma Aldrich S7067-5MG)
- Primocin – 500 mg (InvivoGen Cat. Code: ant-pm-1)
- *R-Spondin (REPROTECH #120-38) optional instead of RCM (100 ng/ml)*
- *Noggin (REPROTECH #250-38) optional instead of NCM (1000 ng/ml)*
- *Wnt3a – if required 100 ng/ml*
- m-EGF (REPROTECH #315-09)
- Real Time – Glo™ MT Cell Viability Assay (Promega #G9711)
- Gastrin (R&D Systems Europe #3006/1)
- Prostaglandin E2 (Sigma-Aldrich #P0409)
- STEMPRO hESC (Fisher A1000701)

1.4

Buffer and Solutions

Buffer preparation for organoids see ADSI – CCO – SOP-1

2.

Reagent preparation

2.1

Media preparation for intestinal organoids

Medium preparation for organoids see SOP-ADSI-1.0.

2.2

Preparation of cells

One well of organoids from a 24 well plate is needed per drug testing condition.

2.2.1

Harvest

- Remove medium from the wells (check if the matrigel is not loose!)
- Add 500 µl cold GF- to every well
- Detach matrigel drop with 1000 µl tip
- Transfer to 15 ml tube containing (keep on ice)
- Wash wells with 500 µl cold GF-
- Add GF- up to 10-14 ml
- Centrifuge for 3 min 800 rpm at 4°C
- Remove supernatant (matrigel pellet with the organoids should be visible by eye).

- If the matrigel pellet is very big, remove the supernatant and resuspend the pellet with 3ml cold GF-. Pipette up and down. Add 7ml of ice cold GF- and centrifuge again for 3 minutes at 1200 rpm at 4°C.


2.2.2 Trypsinisation

- Add 200 to 500 µl 0.25% Trypsin-EDTA (depending on the size of the pellet)
- Pipette up and down
- Incubate 7 min at 37°C (water bath)
- Pipette up and down with a 100 µl pipette tip (if cells will not get loose either extend the incubation time or attach a 200µl pipette tip to the 1000µl pipette tip)
- Check repeatedly under the microscope (in the tube).
- Repeat incubation cycles until the majority of cells are single cells
- If single-cells still do not appear after 3 times of checking under the microscope, add 200 µl 0.25% Trypsin-EDTA and extend incubation cycles at 37°

2.2.2 Seeding

- When single cell suspension is obtained, add 10 ml GF-
- Centrifuge for 5 min, 1500 rpm at 4°C
- Remove the supernatant and resuspend in 1ml of GF- and count the cells.
- Add 10 ml GF-
- Centrifuge for 5 min, 1500 rpm at 4°C
- Resuspend the cells in GF- at 1.5×10^6 cells/ml in the required amount GF- (example: 3,3 µl GF- and 6,7 µl matrigel per well in 96 well plate) resulting in 5000 cells/ well.
- Seed 10 µl droplets per well. Pipette fast to avoid cooling of the pre-warmed plate. (See **seeding scheme**)
- Incubate for 10-15 min at 37°C (5% CO₂)
- Carefully add 150 µl pre-warmed medium (M2 - M6) per well.
- Incubate at 37°C (5% CO₂) for 3 days. (See **time schedule**).
- After 3 days, organoids are ready for the drug screening (Real Time-Glo Viability assay).

Figure 1 Seeding scheme for WT and tumor-organoids on a 96 well plate. Assay to investigate the viability of wild-type and tumor organoids upon drug treatment. Seed 10 μ l droplets of matrigel mixed with GF- and 5×10^3 cells per wells. WT: wild type organoids, TC1: Tumoroids clone 1, TC2: Tumoroids clone 2. The wells F2 to G10 – with matrigel only. Empty wells E2-11 prevent from crosstalk between samples and background controls. Outer wells F2 and G10 (wells with matrigel only) are filled with 150 μ l PBS to prevent evaporation.



Time schedule

Weekday	Day	Media	RTG	Images
Friday	d0	seed in respective EM		
Saturday	d1			
Sunday	d2			
Monday	d3	drug/ vehicle + RTG	RTG timepoint 1h, 2h	+
Tuesday	d4		RTG timepoint 24h	+
Wednesday	d5		RTG timepoint 48h	+
Thursday	d6		RTG timepoint 72h	+

3. Procedures for viability assay conduction

This SOP shows a simple example with 2 different drugs + 1 control (vehicle) tested at 3 different batches of organoid cultures (1x healthy WT organoids + 2 x Tumor organoids in different media –from the same patient). The experimental setup allows the comparison of the anti-tumor activity of a drug on two different organoid clones compared to organoids from healthy tissue.

3.1 Preparatory work

- Put the Real Time-Glo (RTG) Kit at RT, vortex and centrifuge down.
- Prepare 1,25 x of the respective media (#2, #3, #4, #5, or #6). The medium is stable up to 2 weeks at 4°C.
- Table 1 shows media preparation for 4,4 ml 1,25x medium.

- A 5x stock solution of the RTG enzyme and substrate is prepared in GF-. This is essential to guarantee an equally amounts of Real Time-Glo throughout the entire assay.
- The 5x RTG in GF- is than used to prepare the 5x drugs/ vehicle stock solutions (*The 5x solution mixed with the 1,25x medium results in a 1x medium in the well*). The amount of stock solution required depends on how many wells per 96 well plate is to be tested (*always work in triplicates + duplicates for the blank control*). In total a 96 well contains 80 µl of the 1,25 x media solution + 20 µl of the 5 x stock solution.

Media		Expansion Media		M2•Y	M3•Y	M4•Y	M5•Y	M6•Y = WT•Y
Compound	stock conc.	final conc.	dilution	µl	µl	µl	µl	µl
Final volume				5500.00	5500.00	5500.00	5500.00	5500.00
GF-				3379.20	3376.82	3377.37	3374.98	624.43
WCM	1x	50%	2 x	o	o	o	o	2750.00
RCM	10x	10%	10 x	550.00	550.00	550.00	550.00	550.00
NCM	20x	5%	20 x	275.00	275.00	275.00	275.00	275.00
mEGF	10000x	50 ng/ml	10000 x	o	0.55	o	0.55	0.55
Nicotinamide	1M	10 mM	100 x	55.00	55.00	55.00	55.00	55.00
N-acetyl ^l	500mM	1.25 mM	400 x	13.75	13.75	13.75	13.75	13.75
B27	50x	1 x	50 x	110.00	110.00	110.00	110.00	110.00
Y-27632 (ROCK-inh)	10mM	10 µM	1000 x	5.50	5.50	5.50	5.50	5.50
A83-01 (TGFβ-inh)	1.5mM	500 nM	3000 x	o	o	1.83	1.83	1.83
SB202190 (P38 inh)	30mM	10 µM	3000 x	o	1.83	o	1.83	1.83
Primocin	50mg/ml	100 µg/ml	500 x	11.00	11.00	11.00	11.00	11.00
Gastrin	100µM	10 nM	10000 x	0.55	0.55	0.55	0.55	0.55
PGE2	100µM	10 nM	10000 x					0.55
Volume without GF-				1020.80	1023.18	1022.63	1025.02	3775.57
1,25x stock medium				4400.00	4400.00	4400.00	4400.00	4400.00

^lN-acetyl needs to be added last because of it's low pH!

3.2

Drug treatment

- Remove medium from the wells of 96 well plate (do not touch the matrigel drop. Use P1000+P10 tip).
- Carefully add 80 μ l of pre-warmed 1.25x media (Medium #2, #3, #4, #5, #6)
- Add 20 μ l of pre-warmed 5x RTG in GF- containing drugs or vehicle, according to the pipetting scheme.
- Incubate at 37°C (5% CO₂).
- Perform luminescence measurement at time points: 1h, 4h, 24h, 48h, and 72h.
- Take pictures of every well with organoids at time points: 1h, 24h, 48, 72h, (optionally also 96h, 120h, and 144h).
- If you plan to image longer than 72h, than refresh medium after 72 hours to medium without enzyme and substrate, but with the drug.

3.3

Measurement

- Switch on the luminometer (Mitras) at least 1 hour before the measurement.
- Turn on the PC and open the program MikroWin 2000.
- Open program to measure luminescence for 1 and 0.5 s according to your instrument setting. Make sure that the luminometer is at 37°C
- Use a filter-set allowing luminescence measurement
- The measurement is very temperature sensitive. Keep the plates in a styrofoam box for the transfer from the incubator to the luminometer. Try to be as fast as possible to prevent the plate from cooling down.
- Export the data as EXCEL file for further calculations.

4.

References

Broutier L, Andersson-Rolf A, Hindley C J, Boj S F, Clevers H, Koo B K, Huch M. Culture and establishment of self-renewing human and mouse adult liver and pancreas 3D organoids and their genetic manipulation. Nat Protoc. 2016 Sep;11(9):1724-43. doi: 10.1038/nprot.2016.097. Epub 2016 Aug 25.

5.

Applicable references

ADSI – CCO - SOP-1.0, ADSI – CCO - SOP-2.0, ADSI – CCO - SOP-3.0